

Steps for Achieving the Ideal Intraoral Scan:

- Always aim for full arch scans and avoid half arch scans.
- Capture 3-5mm of the patient's buccal and lingual/palatal tissues, including the area behind the last distal tooth.
- If the margin line needs to be scanned, ensure proper gingival retraction, especially with sub-gingival margins. We recommend a two-cord retraction technique for the best visibility.
- Prevent distortion by thoroughly drying the teeth using compressed air, saliva ejectors, and cotton rolls. Keep the interproximal areas dry as well.
- Review the final scan before uploading to check for any missing details or irregularities.

